# **DECATHLON**

## **TW500**

RACKET SPORTS WHEELCHAIR

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Please read this manual carefully before adjusting or using this product.

If necessary, contact Decathlon or the manufacturing company at the address provided below.

Should you suffer from visual impairment, please download the manual in PDF format from our website:

[www.decathlon.com](http://www.decathlon.com)

Wheelchair designed and manufactured by:

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If necessary, and for any further details regarding the safety of the wheelchair, repairs and ordinary or extraordinary maintenance, please contact:

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All rights reserved

# DECATHLON

<b>1. Introduction</b>	<b>4</b>
TW500 Description	4
TW500 Features	5
<b>2. Warning</b>	<b>5</b>
<b>3. Supplied base configuration</b>	<b>6</b>
First use	6
Base configuration	7
<b>4. Adjustments</b>	<b>7</b>
Adjusting the seat depth	7
Adjusting imbalances (COG)	10
Adjusting the height of the fore-aft seat in relation to the ground	10
Adjusting the fore-aft seat height in relation to the ground	12
Adjusting the height of the backrest	12
Adjusting the height of the clothing protection rims	13
Adjusting the backrest incline	13
Adjusting the tautness of the backrest canvas	14
Adjusting the position of the footrest	14
Adjusting the anti-tilt wheel	15
Adjusting the thigh attachment strap	15
Assembling the pelvic strap	16
Adjusting and replacing the hand grip	16
<b>5. Tips and recommendations for use</b>	<b>17</b>
<b>6. Maintenance, inspections and checks</b>	<b>17</b>
Replacing the tyre and air chamber	18
Replacing the castors and bearings	18
Checking the quick removal devices	19
<b>7. Cleaning and disinfection instructions</b>	<b>20</b>
<b>8. Packaging, transport and delivery</b>	<b>20</b>
<b>9. Material identification</b>	<b>21</b>

# 1. Introduction

TW500 is a wheelchair designed to meet the needs of adult and teen tennis and racket sport players. Made with a lightweight aluminium alloy, it offers a wide range of adjustments, meaning it can be configured in line with the athletes' needs and the requirements of their various pathologies. Plus, it means it can accompany them as they develop and progress in their sport, from leisure to competitions. The wheelchair complies with the following standards:


ISO 7176-8 (test model 120 kg)
ISO 7176-16

You must NOT use wheelchairs and their elements incorrectly or for uses other than those indicated in this manual.

## 1.1 TW500 Description

The TW500 tennis wheelchair is characterised by its frame which divides into two sections: A lower frame or base, and an upper frame or seat. The base determines the geometry of the wheelchair's footprint, the seat and the attachment of the drive wheels, the small directing front wheels and the rear anti-tip wheels and includes the footrest and its controls.

- 1. Backrest
- 2. Seat
- 3. Side guards
- 4. Hand grip
- 5. Footrest
- 6. Front guard
- 7. Anti-tip wheels
- RW Rear wheels
- FC Front wheels

## 1.2 TW500 Features

- Aluminium alloy frame.
- D 72 mm height and depth-adjustable anti-tip wheel.
- Predefined camber: 20°
- Drive wheels: 25"
- High pressure tire with 25"x1 (25-559 mm) diameter
- Aluminium push ring
- Rotating front wheels 72 mm diameter, 23.5 mm width
- Selectable seat width: 340, 360, 380, 400 mm
- Seat depth: 340, 360, 380, 400, 420 mm **adjustable**
- Axle/backrest distance: from 60 to 210 mm **adjustable** in 10 mm steps
- Backrest height: 250, 270, 290, 310, 330, 350, 370, 390 mm **adjustable**
- **Backrest incline (in relation to ground): +5° -5°**
- Seat front height in relation to ground: **adjustable** from 400 to 540 mm
- Seat rear height in relation to ground: **adjustable** from 340 to 530 mm in 10 mm increments
- Seat inclination:min -8°, max 34°
- Leg incline in relation to seat: min 50°, max 130°
- Footrest Height, Incline and Distance in relation to the axle: **adjustable**

- Foot strap: **standard**
- Polyurethane cushion 75 sh., H 50 mm (prof. Std. 400 mm) **not anti-sores**.
- Thigh strap H 80 mm (size M 560+ 1100 mm, size L 680+ 1200 mm)
- Shin strap: H 50 mm
- Frame protection: Padding on front and on rear bar
- **Volume:**
  - Total width: min 790 mm, max 930 mm [700 mm]\*\*
  - Total length: min 840 mm, max 890 mm [1200 mm]\*\*
  - Total height: min 570 mm, max 960 mm [1200 mm]\*\*
  - PIVOT width: min 1580 mm, max 1860 mm [1300 mm] \*\*
  - Hand rim: min 790 mm, max 930 mm [1000 mm] \*\*

\*\* Note: Some volume measurements may be greater than the measurements set out in the EN ISO 12183 regulation [xxx mm]. In some cases, it may be difficult or impossible to use emergency exits provided.
- **Weight:** min 15 kg, max 15.9 kg
- **Maximum load:** 120 kg

## 2. Warning

The user is reminded that they should carry out a self-assessment of the critical points associated with their specific pathology, and with their particular postural requirements, before proceeding with any adjustments or use of the wheelchair. **OFFCARR does not accept any responsibility for any damages caused by incorrect seating in the wheelchair, or by features of sitting in the wheelchair which are incompatible with the user's pathology.**

-  The cushion provided **IS NOT ANTI-SORE**, it is made with expanded polyurethane with a density of 75 shore with a Nylon fabric cover.
-  Avoid contact of the wheelchair with water. Undesirable rust may occur on some metal parts and/or loss of the safety features of the materials concerned.
-  The wheelchair must only be used for the sport it has been designed for, or for compatible disciplines, and not for transporting objects in general, and it must not be used in other contexts. This could be dangerous insofar as some of the usual safety elements for the wheelchair in motion are absent (e.g. Parking brakes, etc.).
-  Allergies caused by contact with parts of the wheelchair should be reported.
-  Keep the wheelchair away from heat sources in order to avoid any of the components setting alight. The coverings meet the criteria for the standard ISO 7176 16:2012.
-  Do not put your fingers into the spikes when the wheelchair is in motion.
-  Inflate the tyres to a pressure within the range indicated inside the tyre. If transporting by plane, it is suggested that you deflate the tyres to avoid excess pressure.
-  Always check that the quick wheel removal devices are working properly, especially if you tend to take them off and put them back on frequently.
-  The wheelchair must be regularly maintained, so that it is kept in good, working condition but especially so that the safety conditions are upheld.
-  Inefficient maintenance and the improper use of the device may cause damages and injury to users.
-  Alterations may invalidate the minimum safety conditions.

The registration of the wheelchair and the address of the manufacturer can be found on a small sticker located underneath the wheelchair's lower frame. This sticker must never be removed, at risk of invalidating the product's guarantee.

**The indicated life span of the wheelchair is 2 years under normal usage conditions, subject to it always being used by the same person and undergoing regular maintenance.**


### 3. Supplied base configuration

The TW500 tennis and racket sports wheelchair comes ready-to-use with a standard configuration of the variable parameters.

#### 3.1 First use

The wheelchair is shipped with the rear wheels removed. To position the front wheels, hold down the axle button to add the wheel onto the axle and then push it into the special guide of the wheelchair's frame, then release the button. It is advised that you check the pressure of the wheels before use, ensuring they're at 7.6 BAR. The inflation pressure is indicated on the sides of the tyre covers. **Nevertheless, it is recommended to always check that the drive wheels are inserted correctly before using the wheelchair, namely by checking that the button has been released, and by trying to pull the wheel outwards without pressing on the button.**

#### 3.2 Base configuration

The TW500 wheelchair comes with the following adjustment:

Seat width	L	as requested
Drive wheels	W	25"
Seat depth	D	400 mm
Backrest height	HS	290 mm
Backrest incline	IS	91°
Front height of the seat in relation to the ground	FH	500 mm
Rear height of the seat in relation to the ground	RH	400 mm
Centre of gravity	COG	140 mm
Distance of the footrest from the ground	FD	410 mm
Cushion H 50mm		L x 400 mm
Thigh strap H 80mm		L x 3.5



## 4. Adjustments

### 4.1 Adjusting the seat depth


When selecting the appropriate seat depth, it is recommended that you leave a margin of at least 30 or 40 mm between the front rim of the seat and the back of the knees or the calf. During trial phases, and especially during game play, it is essential that you use a cushion; It is reminded that the cushion provided with the wheelchair is not anti-sore and its depth is 400 mm. If there's no sign or reporting of any particular sensitivity in the skin, you can use the cushion provided. Otherwise, it's preferable to use a cushion adapted to the user's specific needs, such as the cushion they use daily, if it's compatible with the dimensions of the sport wheelchair. Should you have any doubts, please consult your doctor or therapist.

#### 4.1.1 Seat cover depth adjustment


The seat is composed of three separate sections, each with a strap to adjust the tautness.


P 38 (long or short seat frame)


P 42 (long seat frame)

P 36 (short seat frame)

P 34 (short seat frame)

Follow the instructions below to adjust the depth of the seat:


- Remove the "ALS" screws which attach the front section of the "AL" seat to the main frame.
- Assess whether, with the new configuration, you also need to remove the central "AK" section, to

avoid overlaps. In this case, remove the "AKS" screws

- Set the front position in the new position, following the diagram provided, and reinsert the "ALS" attachment screws.
- Put back the central section, if you've removed it, still following the diagram provided, and reinsert the "AKS" screws.
- If necessary, adjust the tightness of the wrist straps and the overlap of the extra fabric pieces on the central section.

#### 4.1.2 Standard configuration of the depth of the seat frame


The depth of the standard frame is 420 mm. This configuration allows you to adjust the seat, or base, to fit the cushion between 340 and 420 mm. It should be emphasised that the curve of the frame is located towards the rim of the seat, when it has a depth between 380 and 420 mm, and is further away or longer when the depth of the seat is between 340 and 380 mm. The images below show the standard configuration.

You should note the position of the "L" telescopic joints with the "xx" split located towards the bottom, and the insertion of the "K" extension between the "J" frame and the "L" joints. The "M" screws attach the telescopic joints to the frame and it's important that they're properly tightened, so that they contribute to the stability of the frame. The "N" screws are solely for attaching the extensions, and for preventing them from rotating so that their upper section is free.

#### 4.1.3 Reducing the depth of the seat frame

This option can only be used when the depth of the seat is adjusted to between 340 and 380 mm.


Please note that, in relation to the standard configuration illustration in paragraph 4.1.2, in this configuration with a reduced depth, the "K" 40 mm extensions are eliminated and the telescopic "L" insert is rotated 180° to allow for the seat to be properly attached. Follow the instructions below to reduce the depth of the frame by eliminating these extensions:

- Firstly, you need to remove the "AL" front section and the "AK" central section of the seat by removing the "ALS" and "AKS" screws as indicated in paragraph 4.1.1.
- Unscrew, but don't remove, the "F1" screws which lock the rotations of the front mount in relation to the base frame (see Fig. Par. 4.2).
- Remove the "F2" screws which lock the rotation of the telescopic joints in relation to the front mounts (take care to also remove the threaded pivot inserted between the rotating joints).
- Remove the "N" screws which attach the K extensions to the joints.
- Remove the "M" screws which attach the telescopic "L" joints to the seat.
- Slide the telescopic "L" joints out of the seat and remove the "K" 40 mm extensions.
- Turn the telescopic "L" joints by 180° by turning the "XX" slot upwards at the start of the rod.
- Completely reinsert the telescopic joints over the seat tube "J".
- Put the "M" screws back in and tighten them properly.
- Assemble the seat, as indicated in paragraph 4.1.1.


## 4.2 Adjusting imbalances (COG)

The imbalance or centre of gravity corresponds to the horizontal distance from the backrest to the axle of the drive wheels. This is a key parameter for distributing the load across the rear drive wheels and the front wheels of the wheelchair, significantly influencing the agility and handling of the wheelchair. The imbalance is set 140 mm in advance, but it can be calibrated to another setting in approximate 10 mm increments:

- Remove the padded protection from the front of the frame.
- Unscrew, but don't remove, the "F1" screws which attach the front mounts to the base of the frame.
- Unscrew, but don't remove, the "F2" screws to release the rotation between the front mounts of the seat.
- Unscrew and remove the "C" screws which attach the seat to the rear support bars.
- Move the seat forwards or backwards to find the desired COG.
- Put back and securely attach the "C" screw in the new position.
- If no other changes are required for the configuration, securely attach the "F1" and "F2" screws which were previously unscrewed.


#### 4.3 Adjusting the height of the fore-aft seat in relation to the ground

The height of the fore-aft seat in relation to the ground is set at 500mm. The configuration provided can be adjusted between 440 and 540 mm. The chair also comes with a kit to accommodate front heights of 420 and 400 mm. To adjust the post position between 440 and 540 mm, follow the instructions below:

- Remove the padded protection from the front of the frame.
- Unscrew, but don't remove, the "F1" screws which attach the front mounts to the base of the frame.
- Unscrew, but don't remove, the "F2" screws to release the rotation of the front mounts in relation to the seat.
- Unscrew and remove the "G" screws which attach the front telescopic sections of the frame.
- Select the desired seat height in relation to the ground, measured between the edge of the seat and the ground (FH).
- Put back and securely attach the "G" screws in the new position.
- If no other changes are required for the seat configuration, properly attach the "F1" and "F2" screws which were previously unscrewed.

To adjust the fore-aft seat position of the wheelchair to 420 or 400 mm, use the hardware provided and follow the instructions below:


- Remove the front padded frame protectors.
- Loosen, but do not remove, the "M" screws that secure the frame curve to the seat
- Unscrew, but don't remove, the "F2" screws to release the rotation between the front mounts of the seat.
- Remove the "F1" screws which attach the "H" interior mounts to the base of the frame.
- Remove the "H" mounts, the "F1" screws, the "D1" spacers and the "G" screws.
- Insert the supplied "A1" rods with the inside of the "F" exterior mounts,
- Attach the "A1" rods to the "F" tube using the "F4" screw, selecting the hole which corresponds to the height desired.
- Attach the exterior mounts to the "M1" bracket using the "F3" screw, "D3" disc washers and "D2" brass sleeves as shown in the illustration.

If no further modification of the seat insert is required, tighten the previously loosened "F2" and "M" screws properly.

#### **4.4 Adjusting the fore-aft seat height in relation to the ground**

The height of the fore-aft seat in relation to the ground can be adjusted to between 340 and 530 mm by following the instructions below:

- Unscrew, without removing, the "F1", "C" and "F2" screws to release any locked rotations and movement of the seat mount.
- Unscrew and remove the "E" screws which attach the rear support bars.
- Select the desired height of the seat in relation to the ground, measured from the end of the seat to the ground "RH".
- Put back the "E" screws in the new position and securely tighten them.
- If no other changes are required for the seat configuration, securely attach the "F1" and "F2" screws which were previously unscrewed.


#### **4.5 Adjusting the height of the backrest**

By default, the height of the backrest is set at 290 mm. Nevertheless, it can be adjusted to between 250 and 390 mm by following the indications below:


- Lift and remove the canvas of the backrest, which is attached by the velcro strips on the support straps.
- Unscrew the (F) screws which attach the adjustment extensions for the backrest height.
- Select the new height for the backrest.
- Put back the (F) screws into the new position and securely tighten them.
- If you need to add one or more velcro support straps for the external padding.
- Re-install the backrest canvas, taking care to fold any overhang under the cushion.

NOTE: The wheelchair comes with two sets of backrest extensions. The short set covers a height adjustment to between 250 and 310 mm. The extended set covers a height adjustment to between 330 and 390 mm. The support straps provided cover the whole adjustment range, just like the padded canvas cover.

## 4.6 Adjusting the height of the clothing protection rims

It's recommended that you define the positioning of the protection height in relation to the wheels before continuing with adjusting the incline of the backrest, in order to eliminate any clashes and needing to redo any adjustments.

- Removing the drive wheels.
- Unscrew and remove the "P", "Q", "V" and "Y" screws which attach the protection to the frame and which set the adjustment of the backrest incline.
- Determine whether you need to remove the incline adjustment. You need to remove it if the "T" screw has also been passed through the protection. In this case, also remove the "T" screw, taking care to note the sequence of the washers for when you put it back together.
- Select the preferred position by using the "P" and "Q" slots nearest to the relevant holes.
- Reinsert all the "P", "Q", "V" and "Y" screws that were previously removed and securely tighten them. If the supporting "T" screw has been removed, it must be reinserted first.


## 4.7 Adjusting the backrest incline

Before adjusting the backrest incline, you should make sure that the clothing protection rim is positioned at the desired height. If necessary, first adjust the height of the rim, following the instructions provided in paragraph 4.7 and then continue as follows with adjusting the incline of the backrest. Operations should always be carried out symmetrically on the right and left sides, in order to be able to make the required rotations. It is recommended that, when the screws are removed, you take note of the sequence of the washers and thicknesses used, so that you can put them back in the right order.

- Unscrew, without removing, the "P" and "Q" screws to release any locked rotations and movements of


the clothing protection rim.

- Unscrew, without removing, the "T" screws to release the rotation of the backrest.
- Unscrew and remove the "V" screw which attaches the incline adjustment onto the seat frame.
- Select the new incline of the backrest.
- Find the hole from those identified by "V1, V2, V3, V4" which match up with one of the holes on the frame, identified by "U1, U2, U3, U4".
- Once you've found the corresponding hole, insert the "V" screw and properly secure it to set the position of the adjustment.
- Continue by tightening the "F", "P" and "Q" screws which were previously unscrewed.


#### 4.8 Adjusting the tautness of the backrest canvas

- Lift the rear section of the canvas.
- Loosen the straps.
- With the athlete in playing position, pull taut the straps as much as necessary, starting with the top strap.
- Fold the backrest canvas back again.
- Try some game movements and repeat the operation if the position isn't quite right.

#### 4.9 Adjusting the position of the footrest

The height, depth and incline of the footrest can be adjusted.

Footplate height is adjusted by varying the positions of the "F5" screws in relation to the support plates, as follows:

- Unscrew and remove the 6 "F5" screws between the plates and blocks "EE" and "FF".
- Align the holes to obtain the desired position.
- Re-insert and properly tighten the "F5" screws.


The angle can be varied by fitting the "F5" screws at different heights to the front and rear plates. Depth is adjusted by varying the position of the "EE" and "FF" blocks in relation to the frame or the platform in relation to the "CC" blocks. Changing the position of blocks "EE" and "FF" is easier and allows for more position possibilities. To change the position of the "EE" and "FF" blocks:

- Loosen the "F6" screws without removing them
- Choose the desired position
- Tighten the "F6" screws firmly

To change the position of the footplate in relation to the "CC" blocks:

- Unscrew and remove the "F7" screws.
- Align the holes to obtain the desired platform position.
- Re-insert and properly tighten the "F7" screws.

## 4.10 Adjusting the anti-tilt wheel

To adjust the depth of the wheel:

- Unlock the "K" collar by pulling the lever
- Increase or decrease the depth of the wheel to suit your requirements
- Tighten the "K" collar
- If the collar is too tight, turn the handle counterclockwise several times, if it is too loose, turn it clockwise until you get the right tension

To adjust the height of the wheel:


- Loosen but do not remove the "AH" screw which determines the anti-tilt wheel's height relative to the ground.
- Raise or lower the wheel as desired.
- Tighten the "AH" screw firmly.

## 4.11 Adjusting the thigh attachment strap

### Preparation

The thigh attachment strap is made up of:

- A padded strap (A) with 2 buckles (B),
- A mobile buckle (C),
- An attachment strap (D).


- Prepare the strap by bringing together both sections "A" and "D" to obtain the desired length ("consider the reference length to be  $W \times 3.25$ ,  $W$  being the width of the wheelchair)
- Insert the strap between the frame and the rims, passing it under the seat.

#### Adjustments


The adjustment of the thigh attachment strap must be made with the athlete sitting in the wheelchair.

- Position strap (A), making sure that the padding is on top of the thighs.
- Insert the attachment strap (D) into one of the two upper buckles (B), taking into account the most comfortable position for the user.
- If the length isn't suitable, it can be adjusted by adapting the connection of the two straps with velcro "A" and "D" which can be found underneath the seat.


#### 4.12 Assembling the pelvic strap

The rim of the TW500 is composed of two elements, including one with slots for potentially attaching a pelvic belt.


#### 4.13 Adjusting and replacing the hand grip

The hand grip can be assembled into three different positions to create a larger or smaller distance "A" from the wheel hoop. To access the hand grip attachment screws, you need to remove the cover, the air chamber [paragraph 6.1] and the protection flap (G).

To remove the hand grip from the wheel, you need to completely remove all six attachment screws (F). Once the hand grip has been removed, you can reattach it by selecting one of the three possible positions on the attachment tabs (see illustration).

Once the hand grip has been re-attached, proceed with the repositioning of the flap (G), making sure to match up the hole for the valve with the corresponding hole in the hoop. We recommend centering the flap (G) over the hoop to ensure maximum protection for the air chamber.

## 5. Tips and recommendations for use

The greater the imbalance, namely the wheelbase between the wheel hub and the backrest axle, the easier the wheelchair is to handle and the better its response to changes in direction. The height of the anti-tip castors in relation to the ground should be adapted to the athlete's game habits. It is recommended that you avoid the possibility of continual contact with the anti-tip castors, insofar as they slow down the wheelchair's movement and interfere with its direction. If needed, the padding can be washed at low temperatures. Nevertheless, avoid immersing the other, mechanical parts of the wheelchair in water. Keep the wheel castors clean, particularly the front and anti-tip castors, and the rotation axles of their forks, so that the wheelchair always functions as efficiently and as effectively as possible.

## 6. Maintenance, inspections and checks

### ATTENTION!

- Always check that the quick removal device for the drive wheels is properly locked.
- The wheelchair should be regularly checked to ensure it remains efficient and safe.
- Insufficient maintenance and improper use of the wheelchair may cause damages and injury to people or property.
- There is no obvious risk of personal injury during use of the wheelchair insofar as it has been designed to reduce the risks of accidents in use or in preparation for use. Nevertheless, as it is a specific product for use in certain sports, it may become hazardous if used in other contexts.
- Should you need to replace any elements, it is recommended that you only use original spare parts or spare parts that have been authorised by the manufacturer.
- Should you encounter any difficulties in finding spare parts on the market, please contact DECATHLON.

Each week:

- Check the tyre pressures. On the side of each wheel cover, you'll find the recommended tyre pressure. A deflated tyre compromises the efficiency of the brakes and the fluidity of the wheelchair's movements.
- Check the efficiency of the quick removal devices [paragraph 6.3], clean them, and if necessary, oil the axles and dowels with a little grease to ensure that the wheels can be simply removed and replaced.

Each quarter:

- Check for wear on the front wheels, which may alter the rear seat of the wheelchair if significantly worn. In this case, proceed with replacing the wheels [par. 6.2].
- Check the bearings on the front and rear wheels, and on the fork rotation. In this case, replace the bearings [par. 6.2].
- Grease the mobile parts, such as the quick removal axles of the wheels and the guides for the wheels. We recommend using silicon oil as a grease, as it's effective and non-staining.

## 6.1 Replacing the tyre and air chamber

### Disassembling the tyre and air chamber

- Deflate the tyre
- Install a special lever between the circle and the shoulder of the tyre, so that when the lever is tilted, the side of the wheel comes out (fig. 1 and 2)
- Install another lever 100 mm from the previous point and repeat the operation (fig. 3)
- Slide the two levers around the whole circle to release the entire tyre.
- Remove the air chamber, starting with the side opposite the inflation valve (fig. 4)
- Once the air chamber has been removed, you can easily remove the tyre to replace either or both of them.

### Assembling the air chamber and tyre

- Remove the air chamber (deflated) in the tyre (fig. 5).
- Insert the air chamber inflation valve into the hole intended for this purpose, located on the circle.
- Around the whole circle, install the tyre, taking note of the direction of the tyre, which indicates whether it's the right or left drive wheel.
- Install the opposite shoulder over the tyre, starting from the point where the valve is, moving round in both directions.
- Insert the final part of the shoulder with the help of the levers intended for this, reversing the instructions shown in Fig. 3, Fig. 2, and Fig. 1.
- Inflate the tyre to the pressure indicated on its side.

## 6.2 Replacing the castors and bearings

### Wheel replacement (Fig. 6)

- Unscrew and remove the "4" screws.
- Remove axle "3" from its seat.
- You can freely remove the wheel.


FIG. 6


FIG. 7

### Replacing the wheel castors (Fig. 7)

- The front wheel castors (608 zz) are mounted with just pressure. To remove them, you just need to push them out from the inside.
- In the illustration, you can see a strut “2” between the two bearings. When reassembling, you must reinsert this strut to set the correct distance between the two bearings.

### Replacing the bearings of the drive wheels (fig. 8)

- Disassembly
  - The drive wheel bearings (6001/12,7 2RS) are mounted by applying pressure. To remove them, you simply need to push them from the inside, using a pin drift and a hammer.
- Assembly
  - Position a new bearing on the inner side of the hub, taking care to insert it evenly (it must be co-axial with the hub, otherwise it won't go in).
  - Before positioning the two bearings, on the opposite side, you need to insert the quick removal axle into the bearing which has already been inserted, so that the axle itself acts as a guide for the two bearings. This way, you make sure that the hub and three bearings are co-axial.


FIG. 8

### 6.3 Checking the quick removal devices

The quick removal axles provided are checked and adjusted: Ø 12.7 mm x Usable L 130 mm. Nevertheless, we recommend carrying out this check regularly during their use.

#### Check:

- Make sure that the “X” button has been entirely released (fig. 10), once the wheel has been assembled.

- Pull the wheel towards the outside of the wheelchair (fig. 11), and the wheel shouldn't come out without you pressing the button.


FIG. 9


FIG. 10


FIG. 11

#### **Adjustment:**

If needed, you can adjust the axle to eliminate any potential alterations between the wheel and the frame, or to complete the release of the button when the wheel has been added.

- If, once the wheel has been fitted into the frame, the axle button isn't completely released, you need to extend the usable length of the axle "L" by partially unscrewing the "Z" nut.
- If, once the wheel has been fitted into the frame, there's interference between the frame and wheel, you need to extend the usable length of the axle "L" by partially unscrewing the "Z" nut.

Note: The "Z" nut has a 1 mm thread, and consequently, the one-turn screwing includes extending or reducing it by 1 mm / If adjusting, it is recommended that you make successive adjustments of 1/4 turn at a time.

## **7. Cleaning and disinfection instructions**

**ATTENTION!** Use the appropriate eye/face protection and gloves during cleaning and disinfection operations.

If the wheelchair comes into contact with blood or other bodily fluids, it must first be cleaned and disinfected, following the instructions below.

Note: In most cases, it is more interesting and more efficient to remove the seat and backrest canvas before cleaning and disinfecting..

#### **Frame**

- Wash the wheelchair with a damp cloth, using warm water and neutral soap to remove dirt or stains.
- Wash the wheelchair again with a damp cloth to remove the detergent product.
- Wipe down the wheelchair
- Visually check that the wheelchair is clean
- disinfect the wheelchair with 70-90% alcohol.
- Make sure that the wheelchair is completely dry before proceeding.

#### **Canvases**

##### If the user is the same before and after cleaning:

- Wash, rinse and disinfect the canvases, using the same procedure as used for the frame.
- Make sure that the canvases are completely dry before re-attaching them to the frame.

##### If the user changes after cleaning:

- We recommend replacing the canvas from the seat and backrest with new, fresh canvases.

NOTE: During cleaning operations, the wheelchair must be carefully inspected for ant damages, rust or malfunctions.. If damage or malfunction is reported, the elements must be repaired or replaced.

**ATTENTION!** All waste generated by this process must be discarded in line with the applicable local legislation.

## 8. Packaging, transport and delivery

All wheelchairs are shipped packaged in sealed cardboard boxes, intended to protect them from impacts and dust. The packaging includes the wheelchair in its basic configuration, fitted with wheels and accessories, a service kit for basic maintenance, and this user manual to which the guarantee card is attached. The wheelchair should be transported using sealed transport methods which protect it from atmospheric agents, as reminded in the instructions provided to this effect on the packaging. Upon delivery, check the integrity of the packaging: Open the packaging and check to make sure that the wheelchair hasn't been damaged in transit. Once these initial checks, which are compulsory for the validity of the guarantee, have been carried out and if the wheelchair isn't used soon after, repack it and store it in a dry place, in a storage area at a temperature between -15 °C and 50° C and RH lower than 80%. Do not lean or place objects on the packaging, as indicated via the symbols featured on the packaging.

## 9. Material identification

The wheelchairs are made using the following materials: aluminium, steel, titanium, polyurethane thermosetting resins, and other composite materials. The sorting and disposal of these materials must be carried out in accordance with the regulations currently in force.

**TW500 products comply with Regulation (EU) MDR 2017 745.**

*OFF CARR s.r.l. Reserves the right to make improvements and/or changes to its products without prior notice, in compliance with the appropriateness of the assistance wheelchair, the contractual conditions of guarantee and the availability of spare parts in compliance with the applicable legal provisions.*

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