DECATHLON

FW500

FENCING WHEELCHAIR
FAUTEUIL ROULANT DE ESCRIME
SILLA DE RUEDAS PARA ESGRIMA
ROLLSTUHL FÜR FECHTSPO
CARROZZINA DA SCHERMA
CADEIRA DE RODAS PARA ESGRIMA
WÓZEK INWALIDZKI DLA SZERMIERZY
MAČEVANJE KOLICA
MEČEVANJE VOZIČEK

INSTRUCTION FOR USE - MODE D'EMPLOI - MODO DE EMPLEO - GEBRAUCHSANWEISUNG - GUIDA PER L'UTENTE - MODO DE UTILIZAÇÃO - INSTRUKCJA OBSŁUGI - UPUTE ZA UPORABU - NAVODILA ZA UPORABO

CE

DECATHLON

ENGLISH	2
FRANÇAIS	17
ESPAÑOL	32
DEUTSCH	47
ITALIANO	63
PORTUGUÊS	78
POLSKIE	93
HRVATSKI	108
SLOVENŠČINA	123

Please read this manual carefully before adjusting or using this product.

If necessary, contact Decathlon or the manufacturing company at the address provided below.

Should you suffer from visual impairment, please download the manual in PDF format from our website:

www.decathlon.com

Wheelchair designed and manufactured by:

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Tel. +39 049 9325733 Fax +39 049 9325734

E-mail: offcarr@offcarr.com <http://www.offcarr.com>

If necessary, and for any further details regarding the safety of the wheelchair, repairs and ordinary or extraordinary maintenance, please contact:

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All rights reserved

DECATHLON

1. Introduction	4
FW500 Description	4
FW500 Features	4
2. Warning	5
3. Supplied base configuration	6
First use	6
Base configuration	6
4. Adjustments	7
Adjusting the height of the backrest	7
Adjusting the seat depth	7
Adjusting the tautness of the backrest canvas	8
Adjusting the distance between the footplate and the seat	8
Adjusting the height and position of the handle	8
Adjusting the thigh attachment strap	9
Assembling the pelvic strap	10
Adjusting and replacing the hand grip	10
Anti-tilt (optional)	10
5. Tips and recommendations for use	11
6. Maintenance, inspections and checks	11
Replacing the tyre and air chamber	12
Replacing the castors and bearings	12
Checking the quick removal devices	13
7. Cleaning and disinfection instructions	14
8. Packaging, transport and delivery	15
9. Material identification	15

1. Introduction

The FW500 is a wheelchair designed to meet the requirements of adults and young people learning fencing on indoor pistes. Made with a lightweight aluminium alloy, it offers a wide range of adjustments allowing it to respond to the athletes' needs and the requirements of their various pathologies. Additionally, it means it can accompany them as they develop and progress in their sport, from leisure to competition. The wheelchair complies with the following standards:

ISO 7176-8 (test model 120 kg)
ISO 7176-16

You must NOT use wheelchairs and their elements incorrectly or for uses other than those indicated in this manual.

1.1 FW500 Description

The FW500 fencing wheelchair features a fixed frame with a large convenient handle on the side of the non-sword arm as well as parking brakes. The height and position (right or left) of the handle, the height of the backrest and the distance between the seat and the footrest are adjustable.

1.2 FW500 Features

- Aluminium alloy frame
- Rim height of sword arm: 80 mm
- Rim height of non-sword arm: 200 mm
- Handle height: 220 mm from the seat (adjustable from 220 to 250 mm)
- Handle depth: 120 mm
- Anti-displacement front bows on the fencing piste
- Camber: 2°
- Rear wheels: 25"
- Aluminium push ramp
- Diameter of road rings for 25" (530 mm)
- Front swivel wheels: D 125 mm
- Selectable seat width: 400 x 460 mm
- Seat depth: 400 mm (adjustable from 400 to 420 mm)
- Distance hub/backrest: 90 mm
- Backrest height: 250, 270, 290, 310, 330, 350, 370, 390 mm **adjustable**
- Backrest angle relative to the ground: 91°
- Seat front height in relation to ground: 530 mm
- Rear seat height from the ground: 530 mm
- Seat angle: 1° fixed
- Leg angle relative to the seat: 83° fixed
- Aluminium backrest
- Polyurethane cushion 75 sh., H 100 mm (prof.. 400 mm) **not anti-sores.**
- Thigh attachment strap: length 1800 mm (L 460 mm) or 1600 mm (L 400 mm) e h 80 mm
- Shin strap: h 50mm
- **Dimensions:**
 - **Total width: min 630 mm, max 690 mm [700 mm]****
 - **Total length: min 860 mm, max 870 mm [1200 mm]****
 - **Total height: min 810, max 940 mm [1200 mm]****
 - **PIVOT width: min 1260 mm, max 1380 mm [1300 mm]****
 - **Turning diameter: min 630 mm, max 690 mm [1000 mm]****
- **Total weight: min 14.3 max 15 kg**

** Note: Some dimension measurements may be greater than the measurements set out in the EN ISO 12183 regulation [xxx mm]. In some cases, it may be difficult or impossible to use emergency exits provided.

- **Maximum user weight: 120 kg**

2. Warning

The user is reminded that they should carry out a self-assessment of the critical points associated with their specific pathology, and with their particular postural requirements, before proceeding with any adjustments or use of the wheelchair.

OFFCARR does not accept any responsibility for any damages caused by incorrect seating in the wheelchair, or by features of sitting in the wheelchair which are incompatible with the user's pathology.

 **The cushion provided IS NOT ANTI-SORE, it is made with expanded polyurethane with a density of 75 shore with a Nylon fabric cover.**

 Avoid contact of the wheelchair with water. Undesirable rust may occur on some metal parts and/or loss of the safety features of the materials concerned.

 The wheelchair must only be used for the sport it has been designed for, or for compatible disciplines, and not for transporting objects in general, and it must not be used in other contexts. This could be dangerous insofar as some of the usual safety elements for the wheelchair in motion are absent (e.g. Parking brakes, etc.).

 Allergies caused by contact with parts of the wheelchair should be reported.

 Keep the wheelchair away from heat sources in order to avoid any of the components setting alight. The coverings meet the criteria for the standard ISO 7176 16:2012.

 Do not put your fingers into the spikes when the wheelchair is in motion.

 Inflate the tyres to a pressure within the range indicated inside the tyre. If transporting by plane, it is suggested that you deflate the tyres to avoid excess pressure.

 Always check that the quick wheel removal devices are working properly, especially if you tend to take them off and put them back on frequently.

 The wheelchair must be regularly maintained, so that it is kept in good, working condition but especially so that the safety conditions are upheld.

 Inefficient maintenance and the improper use of the device may cause damages and injury to users.

 Alterations may invalidate the minimum safety conditions.

 The registration of the wheelchair and the address of the manufacturer can be found on a small sticker located underneath the wheelchair's lower frame. This sticker must never be removed, at risk of invalidating the product's guarantee.

 **The indicated life span of the wheelchair is 2 years under normal usage conditions, subject to it always being used by the same person and undergoing regular maintenance.**

 The brake is for parking purposes only and must not be used as a service brake to slow the wheelchair when in motion.

⚠ Avoid leaning forward excessively, as there is a risk of overturning if the centre of gravity is shifted too far.

⚠ Always negotiate slopes greater than 6° with an assistant.

⚠ Refrain from raising the wheelchair without assistance.

3. Supplied base configuration

The FW500 fencing wheelchair comes ready-to-use with a standard configuration of the variable parameters.




3.1 First use

The wheelchair is shipped with the rear wheels removed. To position the front wheels, hold down the axle button to add the wheel onto the axle and then push it into the special guide of the wheelchair's frame, then release the button. It is advised that you check the pressure of the wheels before use, ensuring they're at 7.6 BAR. The inflation pressure is indicated on the sides of the tyre covers. **Nevertheless, it is recommended to always check that the drive wheels are inserted correctly before using the wheelchair, namely by checking that the button has been released, and by trying to pull the wheel outwards without pressing on the button.**

3.2 Base configuration


The FW500 wheelchair comes with the following adjustment:

Seat width:	L	as required
Main wheel:	W	25"
Seat depth:	D	400 mm
Backrest height:	HS	250 mm
Handle height:	HM	220 mm (from seat)
Handle depth:		120 mm
Rim height of sword arm:	HSA	80 mm
Rim height of non-sword arm:	HNSA	200 mm
Front swivel wheels:		125 mm rotating
Backrest incline:	IS	90°
Front height of the seat in relation to the ground	FH	530 mm
Rear height of the seat in relation to the ground	RH	530 mm
Centre of gravity:	D	90 mm
Distance of the footrest relative to the seat:	FD	350 mm
Cushion H 100 mm		L x 400 mm
Thigh strap H 80 mm		1600 mm if L = 400 mm 1800 mm if L = 460 mm

4. Adjustments

4.1 Adjusting the height of the backrest

By default, the height of the backrest is set at 290 mm. Nevertheless, it can be adjusted to between 250 and 390 mm by following the indications below:


- Lift and remove the canvas of the backrest, which is attached by the velcro strips on the support straps.
- Unscrew the (F) screws which attach the adjustment extensions for the backrest height.
- Select the new height for the backrest.
- Put back the (F) screws into the new position and securely tighten them.
- If you need to add one or more velcro support straps for the external padding.
- Re-install the backrest canvas, taking care to fold any overhang under the cushion.

NOTE: The wheelchair comes with two sets of backrest extensions. The short set covers a height adjustment to between 250 and 310 mm. The extended set covers a height adjustment to between 330 and 390 mm. The support straps provided cover the whole adjustment range, just like the padded canvas cover.

4.2 Adjusting the seat depth


The FW500 allows the seat depth to be increased from 400 to 420 mm simply by moving the seat plate forward by 2 cm:

- Unscrew and remove the seat screws (S).
- Slide the rigid seat forward 2 cm
- Tighten the screws (S)

4.3 Adjusting the tautness of the backrest canvas


- Lift the rear section of the canvas.
- Loosen the straps.
- With the athlete in playing position, pull taut the straps as much as necessary, starting with the top strap.
- Fold the backrest canvas back again.
- Try some game movements and repeat the operation if the position isn't quite right.

4.4 Adjusting the distance between the footplate and the seat

The distance between the footplate and the seat can be adjusted from 350 mm to 430 mm (in 20 mm increments).

- Unscrew and remove the "P" screws.
- Position the footplate at the desired height, aligning the holes in the two tubes
- Replace the "P" screws and tighten them

4.5 Adjusting the height and position of the handle

The FW500 handle is height adjustable.


- Unscrew and remove the "F" thumb screws from the handle brackets
- Raise or lower the handle, ensuring that the holes are correctly aligned
- Replace and secure the thumb screws properly

The FW500 is a wheelchair for right- or left-handed fencers. It is supplied initially with the handle on the left side (for right-handed fencers).

NOTE: Before changing the position of the handle (from left to right or vice versa) it is necessary to switch the rims. To do this, unscrew the screws that fix them to the frame with a 3mm allen key, change their position and re-attach them properly.

To change the handle side :

- Unscrew and remove the "F" thumb screws that secure the handle
- Completely remove the handle and insert it on the opposite side at the desired height, aligning the holes correctly and symmetrically.
- Replace and secure the "F" thumb screws properly


4.6 Adjusting the thigh attachment strap

Preparation

The thigh attachment strap is made up of:

- A padded strap (A) with 2 buckles (B),
- A mobile buckle (C),
- An attachment strap (D).


- Prepare the strap by bringing together both sections “A” and “D” to obtain the desired length “consider the reference length to be $W \times 3.25$, W being the width of the wheelchair)
- Insert the strap between the frame and the rims, passing it under the seat.

Adjustments


The adjustment of the thigh attachment strap must be made with the athlete sitting in the wheelchair.

- Position strap (A), making sure that the padding is on top of the thighs.
- Insert the attachment strap (D) into one of the two upper buckles (B), taking into account the most comfortable position for the user.
- If the length isn't suitable, it can be adjusted by adapting the connection of the two straps with velcro “A” and “D” which can be found underneath the seat.

4.7 Assembling the pelvic strap

The FW500 is designed to be fitted with a pelvic strap attached directly to the frame. The strap is attached to the back tube of the frame and threaded between the rim and the backrest as shown below.


4.8 Adjusting and replacing the hand grip

The hand grip can be assembled into three different positions to create a larger or smaller distance “A” from the wheel hoop. To access the hand grip attachment screws, you need to remove the cover, the air chamber [paragraph 6.1] and the protection flap (G).


To remove the hand grip from the wheel, you need to completely remove all six attachment screws (F). Once the hand grip has been removed, you can reattach it by selecting one of the three possible positions on the attachment tabs (see illustration).

Once the hand grip has been re-attached, proceed with the repositioning of the flap (G), making sure to match up the hole for the valve with the corresponding hole in the hoop. We recommend centring the flap (G) over the hoop to ensure maximum protection for the air chamber.


4.9 Anti-tilt (optional)

The FW500 can be fitted with a device to prevent the wheelchair from tipping backwards. To activate or deactivate it, pull and turn the tube as shown below.

Assembly

- Attach the vice to the mounting plate (Fig. 1)
- Attach the vice to the axle (Fig. 1/2)
- Attach the second vice to the frame (Fig. 2)
- Make adjustments so that the anti-rollover device is parallel to the ground and re-tighten the screws

5. Tips and recommendations for use

If needed, the padding can be washed at low temperatures. Nevertheless, avoid immersing the other, mechanical parts of the wheelchair in water. Keep the wheel bearings clean, especially the front wheels and axles, so that the unit operates as efficiently as possible at all times.

6. Maintenance, inspections and checks

⚠ ATTENTION!

- Always check that the quick removal device for the drive wheels is properly locked.
- The wheelchair should be regularly checked to ensure it remains efficient and safe.
- Insufficient maintenance and improper use of the wheelchair may cause damages and injury to people or property.
- There is no obvious risk of personal injury during use of the wheelchair insofar as it has been designed to reduce the risks of accidents in use or in preparation for use. Nevertheless, as it is a specific product for use in fencing, it may become hazardous if used in other contexts.
- Should you need to replace any elements, it is recommended that you only use original spare parts or spare parts that have been authorised by the manufacturer.
- Should you encounter any difficulties in finding spare parts on the market, please contact DECATHLON.

Each week:

- Check the tyre pressures. On the side of each wheel cover, you'll find the recommended tyre pressure. A deflated tyre compromises the efficiency of the brakes and the fluidity of the wheelchair's movements.
- Check the efficiency of the quick removal devices [paragraph 6.3], clean them, and if necessary, oil the axles and dowels with a little grease to ensure that the wheels can be simply removed and replaced.
- Check the effectiveness of the brakes.

Each quarter:

- Ensure that the rear wheel bearing screws are securely fastened.
- Check for wear on the front wheels, which may alter the rear seat of the wheelchair if significantly worn. In this case, proceed with replacing the wheels [par. 6.2].
- Check the bearings on the front and rear wheels, and on the fork rotation. In this case, replace the bearings [par. 6.2].
- Grease the mobile parts, such as the quick removal axles of the wheels and the guides for the wheels. We recommend using silicon oil as a grease, as it's effective and non-staining.

6.1 Replacing the tyre and air chamber

Disassembling the tyre and air chamber

- Deflate the tyre
- Install a special lever between the circle and the shoulder of the tyre, so that when the lever is tilted, the side of the wheel comes out (fig. 1 and 2)
- Install another lever 100 mm from the previous point and repeat the operation (fig. 3)
- Slide the two levers around the whole circle to release the entire tyre.
- Remove the air chamber, starting with the side opposite the inflation valve (fig. 4)
- Once the air chamber has been removed, you can easily remove the tyre to replace either or both of them.

Assembling the air chamber and tyre

- Remove the air chamber (deflated) in the tyre (fig. 5).
- Insert the air chamber inflation valve into the whole intended for this purpose, located on the circle.
- Around the whole circle, install the tyre, taking note of the direction of the tyre, which indicates whether it's the right or left drive wheel.
- Install the opposite shoulder over the tyre, starting from the point where the valve is, moving round in both directions.
- Insert the final part of the shoulder with the help of the levers intended for this, reversing the instructions shown in Fig. 3, Fig. 2, and Fig. 1.
- Inflate the tyre to the pressure indicated on its side.

6.2 Replacing the castors and bearings

Wheel replacement (Fig. 6)

- Unscrew and remove the "4" screws.
- Remove axle "3" from its seat.
- You can freely remove the wheel.


FIG. 6


FIG. 7

Replacing the wheel castors (Fig. 7)

- The front wheel castors (608 zz) are mounted with just pressure. To remove them, you just need to push them out from the inside.
- In the illustration, you can see a strut "2" between the two bearings. When reassembling, you must reinsert this strut to set the correct distance between the two bearings.

Replacing the bearings of the drive wheels (fig. 8)

- Disassembly
 - The drive wheel bearings (6001/12,7 2RS) are mounted by applying pressure. To remove them, you simply need to push them from the inside, using a pin drift and a hammer.
- Assembly
 - Position a new bearing on the inner side of the hub, taking care to insert it evenly (it must be co-axial with the hub, otherwise it won't go in).
 - Before positioning the two bearings, on the opposite side, you need to insert the quick removal axle into the bearing which has already been inserted, so that the axle itself acts as a guide for the two bearings. This way, you make sure that the hub and three bearings are co-axial.


FIG. 8

6.3 Checking the quick removal devices

The quick removal axles provided are checked and adjusted: \varnothing 12.7 mm x Usable L 130 mm. Nevertheless, we recommend carrying out this check regularly during their use.

Check:

- Make sure that the "X" button has been entirely released (fig. 10), once the wheel has been assembled.
- Pull the wheel towards the outside of the wheelchair (fig. 11), and the wheel shouldn't come out without you pressing the button.


FIG. 9


FIG. 10


FIG. 11

Adjustment:

If needed, you can adjust the axle to eliminate any potential alterations between the wheel and the frame, or to complete the release of the button when the wheel has been added.

- If, once the wheel has been fitted into the frame, the axle button isn't completely released, you need to extend the usable length of the axle "L" by partially unscrewing the "Z" nut.
- If, once the wheel has been fitted into the frame, there's interference between the frame and wheel, you need to extend the usable length of the axle "L" by partially unscrewing the "Z" nut.

Note: The "Z" nut has a 1 mm thread, and consequently, the one-turn screwing includes extending or reducing it by 1 mm / If adjusting, it is recommended that you make successive adjustments of 1/4 turn at a time.

7. Cleaning and disinfection instructions

⚠ ATTENTION! Use the appropriate eye/face protection and gloves during cleaning and disinfection operations.

If the wheelchair comes into contact with blood or other bodily fluids, it must first be cleaned and disinfected, following the instructions below.

Note: In most cases, it is more interesting and more efficient to remove the seat and backrest canvas before cleaning and disinfecting..

Frame

- Wash the wheelchair with a damp cloth, using warm water and neutral soap to remove dirt or stains.
- Wash the wheelchair again with a damp cloth to remove the detergent product.
- Wipe down the wheelchair
- Visually check that the wheelchair is clean
- disinfect the wheelchair with 70-90% alcohol.
- Make sure that the wheelchair is completely dry before proceeding.

Canvases

If the user is the same before and after cleaning:

- Wash, rinse and disinfect the canvases, using the same procedure as used for the frame.
- Make sure that the canvases are completely dry before re-attaching them to the frame.

If the user changes after cleaning:

- We recommend replacing the canvas from the seat and backrest with new, fresh canvases.

NOTE: During cleaning operations, the wheelchair must be carefully inspected for ant damages, rust or malfunctions.. If damage or malfunction is reported, the elements must be repaired or replaced.

 **ATTENTION!** All waste generated by this process must be discarded in line with the applicable local legislation.

8. Packaging, transport and delivery

All wheelchairs are shipped packaged in sealed cardboard boxes, intended to protect them from impacts and dust. The packaging includes the wheelchair in its basic configuration, fitted with wheels and accessories, a service kit for basic maintenance, and this user manual to which the guarantee card is attached. The wheelchair should be transported using sealed transport methods which protect it from atmospheric agents, as reminded in the instructions provided to this effect on the packaging. Upon delivery, check the integrity of the packaging: Open the packaging and check to make sure that the wheelchair hasn't been damaged in transit. Once these initial checks, which are compulsory for the validity of the guarantee, have been carried out and if the wheelchair isn't used soon after, repackage it and store it in a dry place, in a storage area at a temperature between -15 °C and 50° C and RH lower than 80%. Do not lean or place objects on the packaging, as indicated via the symbols featured on the packaging.

9. Material identification

The wheelchairs are made using the following materials: aluminium, steel, titanium, polyurethane thermosetting resins, and other composite materials. The sorting and disposal of these materials must be carried out in accordance with current regulations.

DECATHLON

FW500 products comply with Regulation (EU) MDR 2017 745.

OFF CARR s.r.l. Reserves the right to make improvements and/or changes to its products without prior notice, in compliance with the appropriateness of the assistance wheelchair, the contractual conditions of guarantee and the availability of spare parts in compliance with the applicable legal provisions.

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